



**COTTONDALE VOLUNTEER FIRE DEPARTMENT
INSPECTION / MAINTENANCE CHECKS**

Medical / O2 Bag – Form INSP-00

INSPECTION BY: _____

DATE: _____

TRUCK: **E-25 / R-25**

Side Pouch

CPR Mask
(4) – Pair of XL – Nitrile gloves
(4) – Simple facemasks
Bottle of alcohol based hand sanitizer

Inside Bag

(2) – Convenience bags
(2) – Trauma shears
(1) – Pair of bandage scissors
Ring cutter
Adult BVM w/ mask (check for operation and cleanliness)
Pediatric BVM w/ mask (check for operation and cleanliness)
Infant BVM w/ mask (check for operation and cleanliness)
Yellow disposable blanket
D-sized Oxygen bottle with at least 1000 psi remaining
 Inspect bottle for dents, cracks or any other damage
 Ensure that O-ring is in place and that regulator does not leak
 Bleed bottle to take pressure off regulator
 Ensure that bottle wrench is attached to regulator
Airway kit
 1 each of OPA sizes (orange, red, yellow, green, black, pink)
 1 each of NPA sizes (22, 24, 26, 28, 30, 32)
 (2) – Petroleum based lubricant
 Suction with canister
 French tip attachment for suction
 Adult bulb syringe
 Infant bulb syringe
Green Pouch (Oxygen)
 (2) – Non-rebreather masks
 (2) – Nasal Canula
Blue Pouch (Assessment)
 (2) – Pen lights (check for operation & brightness)
 Blood pressure cuff
 Stethoscope
Yellow Pouch (Environmental Emergencies)
 (2) – Ice Packs
 (3) – Triangular Bandages
 Bottle of sterile water
Red Pouch (Bleeding)
 (10) – 4 X 4 sterile dressings
 (5) – 3 X 3 sterile dressings
 (5) – 2 X 2 sterile dressings
 (10) – Alcohol prep pads
 (3) – Rolls of Kerlix (cling)
 Large roll of clear waterproof tape

**Please email any changes that need to be made to this form to
info@cottondalefirerescue.org**